

## LUNA™ Cell Counting Slides

L12001 (1 box of 50 slides)  
L12002 (10 boxes of 50 slides)  
L12003 (20 boxes of 50 slides)


## Storage

Room temperature

## Product Description

Materials	PS
Dimensions (W x D x H)	75 mm x 25 mm x 2.3 mm
Number of chamber	2 chambers
Chamber volume	10 $\mu$ L

LUNA™ Cell Counting Slides are individually packaged disposable slides used with the LUNA™ cell counter family specifically for bright field cell counting.




Compatible with:

- LUNA™ Automated Cell Counter
- LUNA-II™ Automated Cell Counter
- LUNA-II YFT™ Automated Yeast Cell Counter
- LUNA-FL™ Dual Fluorescence Cell Counter
- LUNA-STEM™ Automated Fluorescence Cell Counter
- LUNA-FX7™ Automated Cell Counter

## Directions for Use

1. Prepare sample to count with or without staining dye.
2. Tear off the top of the wrapping and take out the slide.
3. Load 10  $\mu$ L of sample to one side chamber.
4. Apply the slide to a compatible LUNA™ counter.
5. There are two-chambers for two tests.



### HEADQUARTERS

FL 3  
28 Simindaero 327beon-gil, Dongan-gu  
Anyang-si, Gyeonggi-do 14055  
South Korea

Tel: +82 (31) 478-4185

### USA

7700 Little River Turnpike STE 207  
Annandale, VA 22003  
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

### EUROPE

1 allée Lavoisier 59650 Villeneuve d'Ascq  
France

Tel: +33 (0)3 74 09 44 35

[www.logosbio.com](http://www.logosbio.com)

LBSM-RD-PI-L2-001 Rev.1

## Product Care

- ! Do not press down a slide with or without wrapping.
- ! When handling, hold the edges of the wrapping, not the slide itself.
- ! Widen the opening and gently remove the slide to prevent possible scratches.
- ! Avoid fingerprints on the chamber area.
- ! Do not load any solvents buffer for counting.

## Disclaimer

This product is for research use only.  
Please consult the material safety data sheet for information regarding hazards and safe handling practices.

Additional information is available on our website [www.logosbio.com](http://www.logosbio.com).  
© Aligned Genetics, Inc. 2021. All rights reserved.

## LUNA™ Cell Counting Slides

L12001 (1 box of 50 slides)  
L12002 (10 boxes of 50 slides)  
L12003 (20 boxes of 50 slides)


## Storage

Room temperature

## Product Description

Materials	PS
Dimensions (W x D x H)	75 mm x 25 mm x 2.3 mm
Number of chamber	2 chambers
Chamber volume	10 $\mu$ L

LUNA™ Cell Counting Slides are individually packaged disposable slides used with the LUNA™ cell counter family specifically for bright field cell counting.




Compatible with:

- LUNA™ Automated Cell Counter
- LUNA-II™ Automated Cell Counter
- LUNA-II YFT™ Automated Yeast Cell Counter
- LUNA-FL™ Dual Fluorescence Cell Counter
- LUNA-STEM™ Automated Fluorescence Cell Counter
- LUNA-FX7™ Automated Cell Counter

## Directions for Use

6. Prepare sample to count with or without staining dye.
7. Tear off the top of the wrapping and take out the slide.
8. Load 10  $\mu$ L of sample to one side chamber.
9. Apply the slide to a compatible LUNA™ counter.
10. There are two-chambers for two tests.



## HEADQUARTERS

FL 3  
28 Simindaero 327beon-gil, Dongan-gu  
Anyang-si, Gyeonggi-do 14055  
South Korea

Tel: +82 (31) 478-4185

## USA

7700 Little River Turnpike STE 207  
Annandale, VA 22003  
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

## EUROPE

1 allée Lavoisier 59650 Villeneuve d'Ascq  
France

Tel: +33 (0)3 74 09 44 35

[www.logosbio.com](http://www.logosbio.com)

LBSM-RD-PI-L2-001 Rev.1

## Product Care

- ! Do not press down a slide with or without wrapping.
- ! When handling, hold the edges of the wrapping, not the slide itself.
- ! Widen the opening and gently remove the slide to prevent possible scratches.
- ! Avoid fingerprints on the chamber area.
- ! Do not load any solvents buffer for counting.

## Disclaimer

This product is for research use only.  
Please consult the material safety data sheet for information regarding hazards and safe handling practices.

Additional information is available on our website [www.logosbio.com](http://www.logosbio.com).  
© Aligned Genetics, Inc. 2021. All rights reserved.